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14.10.04 CB/GP/TM

Operating Manual Mode d'emploi

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Deutsch

Tubing pump
with fix speed drive

MS-CA
Stand-mounted pump
with 2 or 4 channels
6 or 8 rollers

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Contents

Safety precautions	4
Warranty terms	7
Product	8
Package contents	9
Rear panel	10
Mains voltage	10
Changing the fuses	11
Operating panel	
Starting the pump	12
Inserting the tubing	14
Pumping against back-pressure	15
Running in period for tubing	16
Tubing live	16

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Contents

Flow rates	17
Accessories	18
Maintenance	19
Repairs	19
Spare-parts	19
Technical specifications	20

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Please note

We recommend you to read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC SA does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC SA does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® tubing pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- The pump **must not** be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes

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Safety precautions

- The pump must be switched off when cassettes or tubing are inserted or changed.
- The permeability of pumptubing depends on the material used and pressure conditions. Tubing can also become electrostatically charged. Please be aware of possible hazards when laying tubing in explosionproof chambers.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- Tubing pumps consist of revolving parts (rotor of the pumphead). Therefore, the pump must not be operated before the cassettes are fully snapped-in.
- Do not manipulate the pumphead before the pump is switched off and disconnected from the mains.

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Safety precautions

- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.
- Only new fuses, according to the specifications stated on page 10 in this manual, must be used. The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC SA denies any responsibility.

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Warranty

2 years from date of delivery

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

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Warranty Terms

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion. A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality. For pump-heads from manufacturers other than ISMATEC SA the warranty terms of the specific manufacturer are valid. The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded. Shipping costs are charged to the customer.

Our warranty becomes invalid in case of

- improper operation by the user, or if the pump is diverted from its proper use,
- unauthorized modification or misuse by the user or by a third party,
- improper site preparation and maintenance,
- operation outside environmental and electrical specifications for the product,
- use of third-party software, hardware, interfacing or consumables purchased by the user and which do not comply with our specifications,
- damages caused by contamination or leaks due to torn or bursted tubing

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Product

The tubing pump is with 2 or 4 tube channels and 6 or 8 pump rollers available.

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Product

Package contents

- MS-CA pump (as ordered)
including Click'n'go Cassettes
- Power cord
- Operating instruction

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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

^{*)} slow-blow type

⚠ Socket/Power cord

Use exclusively the originally supplied power cord. The socket must be earthed (protective conductor contact).

⚠ Before starting-up

Check if the voltage setting rating on the type plate complies with your local mains voltage.

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
Changing the fuses

→ Switch the pump off,
pull out the mains plug.

1 Pull out the fuse-holder by opening
it with a small screw-driver (size 0).

2 Insert 2 new fuses
230 V_{AC}: 2 x 100 mA (slow-blow)
115 V_{AC}: 2 x 200 mA (slow-blow)

⚠ Use always 2 slow-blow fuses
of the same type complying with
the local mains voltage.

3 Shut the fuse-holder.

⚠ Before starting-up

Check if the voltage setting of the
fuses complies with your local mains
voltage.

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
For applications with a high differential pressure (>1 bar), we recommend the use of cassettes with pressure lever.

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Operating panel

- 1 Mains socket
- 2 Fuse holder
- 3 Switch for Run/Stop and flow direction
Position A = clockwise direction
Position B = stop
Position C = counter-clockwise direction

Starting the pump


- 1 Check the voltage rating on the type plate (on rear panel). If necessary, change the rating as indicated on page 10–11.
- 2 Connect the pump tubing to the system.
- 3 Connect the pump to the mains.

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Starting the pump

4. Depending on the required flow direction, press the upper or lower part of the Run/Stop switch. If the pump is correctly connected to the power supply, the pump-head now starts to rotate in the required direction.

5. For stopping the pump, put the rocker switch back to its neutral middle position "B".


⚠ When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case we recommend you to wet the tubing and to run the pump with the tubing inserted for about 15 to 30 minutes.

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Ersatz-Kassetten MS/CA Click'n'go
Bestell-Nr. IS 3510

Inserting the tubing

- 1 ➔ Switch the pump off
- 2 Remove the cassette by slightly pressing the fixing-tongue and lifting it simultaneously
- 3 Insert the 3-stop color-coded tubing with one collar into the cassette
- 4 Let the tubing hang down (prevent it from being twisted)
- 5 Insert the tubing with the second collar at the other end of the cassette
- 6 Reinsert the cassette into the roller-head
- 7 ➔ Switch the pump on

➔ When the pump is idle,
release all cassettes at the fixing-
tongue (2).

Spare-cassettes MS/CA Click'n'go
Order-No. IS 3510

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Cassettes suppl. MS/CA Click'n'go
No de commande IS 3510

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Pumping against pressure

The MS-CA pump can be used for continuous duty at a differential pressure of max. 1.0 bar (smaller tube sizes and/or lever type cassettes allow higher pressures)

When the pump is idle

When the pump is idle, we recommend you to release the tubing from pressure. Releasing the cassette on the right side is sufficient. This helps to protect the tubing from unnecessary strain and prolongs its service-life.

Syphoning effect

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

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⚠ Please note

Please refer to our detailed tubing documentation.

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Running-in period for tubing

Every new tube requires a running-in period. If constant and reproducible flow rates are required, we recommend you to run new tubing in with water or the medium to be pumped for at least 1 to 3 hours before you start the application.

Tubing life

The service-life of the tubing depends on the application and the tubing material used.

Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated in our ISMATEC® catalog (without obligation!).

In order to improve the lubrication and service-life of the tubing, we recommend users to spray both the tubing and the pump rollers with our silicone oil spray (Order No. SC0179).

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⚠ Flow rates

In the table listed subsequently the values indicated are only approximate and determined as follows:

ml/min, per channel, with water and Tygon® ST tubing, without differential pressure

We will be pleased to help the user to select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you to use tubing from ISMATEC®

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1) Can only be used as single-channel pump

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MS/CA Click'n'go Kasette
MS/CA Click'n'go cassette

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MS/CA pressure lever*

Order No. IS 3629

* The cassettes with pressure lever are available on request. Depending on the application, this type of cassette may provide better results at elevated differential pressure conditions.

Contents of tubing selection guide:

- Material and tubing properties
- Chemical resistance chart
- Ordering information

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Maintenance

Provided the tubing pump MS-CA is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

Repairs

For repairs please send the defective MS-CA to your ISMATEC® representative.

Please give information on defect, date of purchase, serial-no., and model.

Replacement parts

For repairs you intend to carry out by yourself (out of the warranty period), you can ask your ISMATEC® distributor for:

- parts lists
- replacement parts for pump drive
- wiring diagrams

Please give information on defect, date of purchase, serial-no., and model.

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Motor type Synchronous motor
 Fixed speed

Speed ranges:

2- and 4 channels 20 rpm
only 4-channels 40 rpm
only 2-channels 60 rpm

Differential pressure

max. 1.0 bar (14.5 psi)
Depending on tubing material.
Tubing with small i.d.'s and/or
pressure lever cassettes may enable
higher pressures.

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Please note

We also recommend you to observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Technical Specifications

Operating conditions

Temperature +5 to +40°C,
Rel. humidity max. 80%
– not condensing, at normal environmental conditions

Dimensions/Weight

	(DxWxH)
2 channels	125x88x135 mm (4.9x3.5x5.3 inch)
4 channels	145x88x135 mm (5.7x3.5x5.3 inch)

	Weight
2 channels	1.2 kg (2.65 lb)
4 channels	1.3 kg (2.87 lb)

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Ask for the latest
ISMATEC-catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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A Unit of IDEX Corporation
Laboratoriumstechnik GmbH